

# High Vacuum Gas Sorption Analyzers

**Autosorb Series**





1



2



3



4




5

# The Choice Is Clear **Absolutely Autosorb**

With the new Autosorb series, customizable high vacuum gas sorption analyzers designed for the most challenging measurements of BET surface area, active area, and pore size distributions in the nanometer range, we deliver on the most important requirements of material characterization laboratories: ACCURACY, AGILITY, ACCESSIBILITY, ADAPTABILITY, and ASSURANCE.

FIND OUT MORE



[www.anton-paar.com/  
apb-autosorb](http://www.anton-paar.com/apb-autosorb)

1

#### ABSOLUTELY ACCURATE

- Precise manifold temperature control <math><0.05\text{ }^\circ\text{C}</math>
- Exceptionally leak-tight system
- TruZone active coolant level control

2

#### ABSOLUTELY AGILE

- Six best-in-class degassing stations
- Three samples, three gases, three temperatures simultaneously
- 90+ hour cryogen Dewar or 1,100 °C furnace

3

#### ABSOLUTELY ACCESSIBLE

- DoseWizard
- PowderProtect
- Intuitive Kaomi software

4

#### ABSOLUTELY ADAPTABLE

- Three instrument models
- 8+ factory-installed options
- 7+ modular field upgrades

5

#### ABSOLUTELY ASSURED

- Compliant with 20+ ASTM, DIN, and ISO standards
- 3-year warranty
- Anton Paar global support network

# Absolutely Accurate

## Free from Error, Exact

Leveraging over 50 years of experience in gas sorption technology, the Autosorb series has been engineered from the ground up to optimize every step in the measurement process for accuracy and performance.

Autosorb 6100

1

### Precise manifold temperature control

Control the manifold temperature between 35 °C and 50 °C with stability better than 0.05 °C. Obtain accurate results day in and day out, no matter the environmental conditions.

2

### Exceptionally leak-tight system


Vacuum-brazed, stainless steel manifold block construction and pneumatically actuated bellows valves ensure 38 % better vacuum performance than previous-generation instruments and high-accuracy data in the  $10^{-8}$   $p/p_0$  domain with nitrogen.

3


### TruZone active coolant level control

TruZone active coolant level control generates a near-constant temperature profile along the analysis cell ensuring accurate measurement of a wide range of material properties and sample types. The small and constant cold zone also ensures better than 2 % reproducibility of nitrogen BET surface areas even with less than 2 m<sup>2</sup> in the cell.


1



2



3



# Absolutely Agile

## Able to Move Quickly and Easily

Best-in-class analysis and sample preparation features on the Autosorb series keep your research agile, maximizing your laboratory's output of advanced measurements and novel materials.



### Up to three independent analysis stations keep you flexible

Because the analysis stations are separate and have dedicated transducer sets, you can analyze up to three different samples with three different analysis gases simultaneously. You can also pair these with our patented CryoSync accessories so that each station measures at independent analysis temperatures.



### 90+ hour analysis Dewar for long-lasting measurements

The 3-liter Dewar included with Autosorb instruments holds enough cryogen to last for more than 90 hours. You'll never have to come in on the weekend to refill the Dewar again. Connect your Autosorb with alternative Dewars and external temperature control accessories for maximum analysis flexibility.



### Measurement agility realized with quick swap to 1,100 °C furnace

Switch between the cryogenic Dewar and the 1,100 °C furnace included with Autosorb 6200 and 6300 instruments in minutes. Rapidly perform a complete characterization of your catalysts, including BET surface area, pore size distribution, active area, and metal dispersion, quickly and easily.



### Six degassing stations, two independent heating zones, maximum flexibility

Simultaneously prepare two sets of samples at two independent temperatures up to 450 °C. A built-in 2-liter cold trap, turbomolecular pump, and adaptive heating profiles, including automated pressure-controlled heating and test for completion, ensure your samples are properly prepared for analysis.

# Absolutely Accessible

## Easy to Understand, Easy to Use

Take advantage of the new Kaomi software, which gives you streamlined control of your Autosorb – no matter your experience level – while also offering advanced analytical performance and flexibility.



### **DoseWizard system for simplified workflows**

Kaomi for Autosorb makes gas sorption easy for users of all levels. Simply select one of our 45+ preloaded analysis profiles or tell the instrument your desired analysis methods, and Autosorb will do the rest. For expert users, an advanced mode gives full control of analysis procedures.

### **PowderProtect keeps instrument downtime to a minimum**

Our new PowderProtect feature virtually eliminates the risk of sample elutriation (loss of powder into the instrument) without requiring special evacuation parameters for different sample types, preventing contamination of the instrument by inexperienced operators.

### **Sample ID tracking ensures traceability**

The Kaomi for Autosorb sample library stores and tracks all sample IDs, sample weights, cell selections, and degassing conditions so you don't have to.

### **Streamlined user interface keeps you in control**

Monitor the instrument status from the sidebar, log view, or an enlarged instrument schematic view. Ensure the top performance of your instrument by tracking the next service date and performing fully guided maintenance routines. With Kaomi, you're always in control.

### **Dynamic Kaomi React interface makes complex analyses easy**

Analyzing dynamic flow chemisorption measurements has never been easier. The new Kaomi React interface helps you easily find measurement files and intuitively guides you through advanced peak deconvolution and pulse titration analyses.

### **Multilanguage support for operators around the globe**

With support for eight different languages, Kaomi for Autosorb lets you control the instrument and analyze the data in the language in which you are most comfortable.

# Absolutely Adaptable

## Adjusts to New Conditions, Flexible

Customize the Autosorb to meet your current application needs with a choice of three instrument models. As your research changes, modular field upgrades allow your Autosorb to adapt to new applications.

### Autosorb 6100


Our basic high-vacuum analyzer features independent analysis stations and a long-lasting cryogen Dewar perfect for analyses of surface area and pore size. Field upgrades let it adapt to new applications as your needs change.

### Autosorb 6200

Our mid-range analyzer can be customized to perform a wide range of chemisorption and physisorption analyses. Along with the basic features of the Autosorb 6100, this model includes a 1,100 °C furnace, a choice of chemical compatibility, and additional measurement options for advanced catalyst characterization.

### Autosorb 6300

Our most powerful instrument features: full chemisorption and physisorption capabilities ideal for advanced catalyst characterization applications. In addition, PFE elastomers offer the highest chemical resistance of the three instruments and give you access to the widest range of potential measurements in the series.



		Autosorb 6100	Autosorb 6200	Autosorb 6300
		↓	↓	↓
INSTRUMENT FEATURE	APPLICATION			
	Standard resistance (N <sub>2</sub> , Ar, Kr, CO <sub>2</sub> , H <sub>2</sub> , )	✓	✓	✓
Sealing materials / chemical compatibility	CH <sub>4</sub> , C <sub>2</sub> H <sub>6</sub> , C <sub>3</sub> H <sub>8</sub> , and other alkanes, C <sub>6</sub> H <sub>6</sub> , C <sub>8</sub> H <sub>10</sub>	✓	○	✓
	NH <sub>3</sub> , C <sub>5</sub> H <sub>5</sub> N, C <sub>3</sub> H <sub>6</sub> O, and other reactive gases		○	✓
Turbo pump system and low pressure transducers	BET surface area Pore size distribution	✓	✓	✓
Vapor option	Water activity Moisture uptake	○ / U	○ / U	✓
High-temperature (1,100 °C) furnace	Active area		✓	✓
	Dispersion	U		✓
	Crystallite size			✓
Integrated TCD	Reduction/oxidation			
	Acid site strength		○ / U	✓
	Activation energy			
Integrated mass spectrometer	Reactive species identification		○ / U	○ / U

- ✓ Included in base instrument
- Factory-installed option
- U Modular field upgrade

# Absolutely Assured

## Protected, Secure, Confident

With compliance to 20+ ASTM, DIN, and ISO standards, you can be certain that the Autosorb series is the right choice for your characterization needs. What's more, Anton Paar's global network of subsidiaries and distribution partners ensures a qualified expert is always nearby and ready to help.

### 1 Advanced Ceramics

- ASTM C1274 Advanced ceramic specific surface area
- ISO 18757 Specific surface area of ceramic powders

### 2 Carbons


- ASTM D6556 Total and external surface area

### 3 Catalyst and Catalyst Carriers

- ASTM D3663 Surface area
- ASTM D4222 Nitrogen adsorption and desorption isotherms
- ASTM D4365 Micropore volume and zeolite area
- ASTM D4641 Pore size distributions
- ASTM D4780 Low surface area
- ASTM D3908 Hydrogen chemisorption
- ASTM D4824 Catalyst acidity by ammonia chemisorption

### 4 Other Non-Porous and Nanoporous Solids

- ASTM B922 Metal powder specific surface area
- ASTM C1069 Specific surface area of alumina or quartz
- ASTM D1993 Precipitated silica – surface area
- DIN 66134 Pore size distribution and specific surface area
- DIN 66135, 1-4 Particle characterization – micropore analysis
- ISO 9277 Specific surface area of solids
- ISO 15901, 2-3 Pore size distribution and porosity of solid materials



# Reliable. Compliant. Qualified.

Our well-trained and certified technicians are ready to keep your instrument running smoothly.

FIND OUT MORE



[www.anton-paar.com/service](http://www.anton-paar.com/service)



Maximum uptime



Warranty program



Short response times



A global service network

	Autosorb 6100	Autosorb 6200	Autosorb 6300
	↓	↓	↓
<b>MEASUREMENT SPECIFICATIONS</b>			
<b>Measurement principle</b>	- Vacuum Volumetric	- Vacuum Volumetric - Dynamic Flow (Optional)	- Vacuum Volumetric - Dynamic Flow
<b>N<sub>2</sub>, Ar, Kr, CO<sub>2</sub>, H<sub>2</sub>, and other non-corrosive gases</b>	Included		
<b>CH<sub>4</sub>, C<sub>2</sub>H<sub>6</sub>, C<sub>3</sub>H<sub>8</sub>, other alkanes, C<sub>6</sub>H<sub>6</sub>, C<sub>8</sub>H<sub>10</sub></b>	Included	Dependent on configuration	Included
<b>NH<sub>3</sub>, C<sub>5</sub>H<sub>5</sub>N, C<sub>3</sub>H<sub>6</sub>O, and other reactive gases</b>	N/A	Dependent on configuration	Included
<b>Physisorption analysis stations</b>	- Number: 1, 2, or 3 - Independence: Up to 3 gases at 3 analysis temperatures can be used concurrently, 1 analysis gas and temperature per station		
<b>Chemisorption analysis stations</b>	N/A	1	
<b>Independent p<sub>0</sub> station</b>	Yes (dedicated cell and transducer)		
<b>Pressure</b>	- Range: 2x10 <sup>-5</sup> to 1100 Torr (2.6x10 <sup>-8</sup> to 0.997 p/p <sub>0</sub> for N <sub>2</sub> 77K) - Resolution (MP): 2x10 <sup>-5</sup> Torr (2.6x10 <sup>-8</sup> p/p <sub>0</sub> for N <sub>2</sub> 77K) - Resolution (XR): 1x10 <sup>-6</sup> Torr (1.3x10 <sup>-9</sup> p/p <sub>0</sub> for N <sub>2</sub> 77K)		
<b>BET surface area</b>	- Absolute detection limit: 0.1 m <sup>2</sup> (N <sub>2</sub> 77K) - Specific detection limit: 0.01 m <sup>2</sup> /g (N <sub>2</sub> 77K) - Typical reproducibility: 1 % (measured on BAM P115) - Reproducibility limit: 2 % with 2 m <sup>2</sup> in the cell		
<b>Pore size</b>	- Range: 0.35 nm to 500 nm (diameter) - Typical reproducibility: 0.5 % (measured on BAM P115)		
<b>Active area</b>	N/A	- Absolute detection limit: 0.03 m <sup>2</sup> (H <sub>2</sub> on Platinum 313 K) - Specific detection limit: 0.003 m <sup>2</sup> /g (H <sub>2</sub> on Platinum 313 K) - Typical reproducibility: 2 % (measured on 2 % Platinum on Alumina reference sample)	
<b>TruZone</b>	Yes (active coolant level control)		
<b>PowderProtect</b>	Yes (prevents sample elutriation)		
<b>Analysis Dewar</b>	- Duration: 90+ hours with liquid nitrogen - Refill during analysis: Yes - Volume: 3 L		
<b>Analysis furnace</b>	N/A	- Maximum temperature: 1,100 °C - Maximum Ramp Rate: 50 °C per minute - Furnace cooling by built-in fan	
<b>Sample preparation</b>	- Integrated degassing stations with dedicated cold trap: 6 - Temperature control: 2 independent heating zones, ambient to 450 °C - Available methods: flow and vacuum, programmable multi-step degassing profiles, test for completion routines, pressure-controlled heating		

<b>ADDITIONAL SPECIFICATIONS</b>			
<b>Kaomi for Autosorb software</b>	- Control up to 4 instruments from 1 PC - 8 languages: Chinese, English, French, German, Japanese, Korean, Portuguese, Spanish - DoseWizard: 45+ built-in analysis profiles (ASTM, USP, DIN, ISO)		
<b>Vapor source</b>	Optional	Optional	Included
<b>CryoSync control</b>	Optional	Optional	Optional
<b>TCD and loop injector</b>	N/A	Optional	Included
<b>Mass spectrometer</b>	N/A	Optional	Optional

	Autosorb 6100	Autosorb 6200	Autosorb 6300
	↓	↓	↓
<b>TECHNICAL SPECIFICATIONS</b>			
<b>Dimensions (W x D x H)</b>	79.5 cm x 70.1 cm x 107.9 cm (31.3 in x 27.6 in x 42.5 in) dependent on configuration		
<b>Weight</b>	136.4 kg (300 lbs) dependent on configuration		
<b>Operating environment</b>	- Temperature: 15 °C to 35 °C (59 °F to 95 °F) - Humidity: 20 % RH to 80 % RH, non-condensing - Indoor use only		
<b>Gas distribution manifold elastomers</b>	FKM	FKM or EPDM (dependent on configuration)	PFE
<b>Analysis manifold construction</b>	- Vacuum-brazed manifold block in stainless steel		
<b>Analysis manifold temperature</b>	- Range: user adjustable 35 °C to 50 °C - Stability: ±0.05 °C		
<b>Analysis manifold valves</b>	- Type: stainless steel, air-operated bellows valve - Cycle life: 5,000,000 cycles - External leak rate: 5x10 <sup>-12</sup> Pa m <sup>3</sup> /sec - Seat leak rate: 5x10 <sup>-12</sup> Pa m <sup>3</sup> /sec		
<b>Pressure transducer (1550 Torr) - Analysis and degas</b>	- Number: 3 to 5 (dependent on configuration) - Accuracy: ±0.1 % of span typical, ±0.15 % of span maximum - A/D data acquisition: 32 bit		
<b>Pressure transducer (10 Torr) - MP or XR</b>	- Number: 1 to 3 (dependent on configuration) - Accuracy: ±0.15 % reading up to 10 Torr - Resolution: 0.002 % full scale		
<b>Pressure transducer (1 Torr) - MP only</b>	- Number: 1 to 3 (dependent on configuration) - Accuracy: ±0.15 % reading up to 1 Torr - Resolution: 0.002 % full scale		
<b>Pressure transducer (0.1 Torr) - XR only</b>	- Number: 1 to 3 (dependent on instrument configuration) - Accuracy: ±0.15 % reading up to 0.1 Torr - Resolution: 0.001 % full scale		
<b>Pirani vacuum gauge - Degas only</b>	- Number: 1 - Accuracy: 5x10 <sup>-4</sup> to 1x10 <sup>-3</sup> Torr: ±10 % of reading 1x10 <sup>-3</sup> to 100 Torr: ±5 % of reading 100 Torr to atmosphere: ±25 % of reading		
<b>Vacuum system</b>	- Turbomolecular drag pump and dry diaphragm pump - Manufacturer's specification: 5x10 <sup>-10</sup> mbar - Typical pressure at analysis port: 2.67x10 <sup>-5</sup> mbar		
<b>Thermal conductivity detector (TCD)</b>	N/A	Two-filament TCD with rhenium/tungsten filaments Detection limit: 0.5 uL (air injected into helium)	
<b>Included gas ports</b>	7 (5 analysis, 1 helium, 1 degas/backfill)		15 (12 analysis, 1 helium, 1 degas/ backfill, 1 titration)
<b>Optional gas ports</b>	7 additional analysis	7 additional analysis 1 titration gas	N/A
<b>Gas supply</b>	- Purity: 99.999 % - Input pressure: 8 psig to 10 psig (0.55 bar to 0.69 bar)		
<b>Air</b>	Input pressure: 50 psig to 100 psig (3.45 bar to 6.9 bar)		
<b>Electrical</b>	- Supply: 100 to 240 VAC ±10 %, 50 or 60 Hz - Maximum consumption: 1600 W (dependent on configuration)		
<b>PC connection</b>	Ethernet		
<b>RoHS3 compliant</b>	Yes		
<b>CE / UKCA certified</b>	Yes		



