

A Chemist's Guide to Sample Preparation

Basics of Microwave-assisted Acid Digestion in a Nutshell – How to Get Started with Anton Paar Instrumentation



Markus Michaelis, Alexander Städler

A Chemist's Guide to Sample Preparation

Basics of Microwave-assisted
Acid Digestion in a Nutshell –
How to Get Started with
Anton Paar Instrumentation

by
Markus Michaelis,
Alexander Stadler

- Copyright©2021 by Anton Paar GmbH, Austria.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form by any means electronic, mechanical, photocopying or otherwise without first obtaining written permission of the copyright owner.

- Published by **Anton Paar** GmbH. Printed in the EU.

Anton Paar GmbH

Anton-Paar-Str. 20

A-8054 Graz

Austria - Europe

Tel.: +43 316 257-0

Fax: +43 316 257-257

E-Mail: info@anton-paar.com

Web: www.anton-paar.com

2nd Edition, August 2021

1st Edition, January 2021



Contents

1.	Introduction	7
1.1.	About the Book	7
1.2.	Anton Paar and Microwaves	7
1.2.1.	<i>Alternative Sample-Preparation Equipment</i>	9

Part A – How To get Started with Acid Digestion

2.	Basics of Acid Digestion	10
2.1.	Why Acid Digestion?	10
2.2.	Commonly Used Reagents	11
2.2.1.	<i>Specific Aspects of Hydrofluoric Acid</i>	15
2.3.	Types of Transfer into Liquid Form	16
2.3.	Error Sources in Sample Preparation	17
2.4.1.	<i>Contamination</i>	17
2.4.2.	<i>Loss of Analyte</i>	19
2.4.3.	<i>Incomplete Digestion Processes</i>	20
2.5.	Typical Analytical Methods	21
2.5.1.	<i>ICP-MS</i>	21
2.5.2.	<i>ICP-OES</i>	22
2.5.3.	<i>MIP-OES</i>	23
2.5.4.	<i>AAS</i>	24
2.6.	Why Microwave Heating?	25
2.6.1.	<i>Dielectric Heating</i>	26
2.6.2.	<i>Benefits of Microwave Heating Compared to Conventional Heating</i>	28
2.6.3.	<i>Microwave Cavity Variations</i>	29
2.7.	It's All About the Temperature	32
2.8.	Dealing with the Reaction Pressure	34
2.8.1.	<i>Vapor Pressure</i>	34

2.8.2.	<i>Generation of reaction gases</i>	36
2.8.3.	<i>Estimation of Vessel Pressure</i>	38
2.8.4.	<i>Pressure Control</i>	39
2.8.4.1.	Sealed vessels	39
2.8.4.2.	Closed vessels with venting capability	40
2.8.5.	<i>Limitations from Vapor Pressure</i>	41
2.8.5.1.	Sealed vessels	42
2.8.5.2.	Closed vessels with venting capability	42
2.8.6.	<i>Correlation between Sample Weight and Pressure</i>	43
2.9.	Understanding the Use of Sensors	44
2.9.1.	<i>Temperature Sensors</i>	44
2.9.2.	<i>Pressure Control</i>	46
2.10.	Power vs. Temperature Programming	46
2.11.	How to Choose the Right Instrument Configuration	47
2.11.1.	<i>Sealed Vessels vs. Closed Vessels with Venting Capability</i>	49
2.11.2.	<i>Vessel material</i>	50
2.11.3.	<i>Vial Types in Autoclave Instruments</i>	51
3.	Method Development Routines	52
3.1.	Access Application Knowledge	53
3.2.	Working with Generic Methods	56
3.2.1.	<i>Reactive Samples</i>	56
3.2.2.	<i>Inert Samples</i>	60
3.2.3.	<i>Using Standard Methods</i>	62
3.3.1.	<i>Increase Oxidation of Organic Matter</i>	63
3.3.2.	<i>Decrease Reactivity of HNO₃</i>	64
3.3.3.	<i>H₂O₂ Addition – Why and How Much?</i>	64
3.3.4.	<i>Optimizing Parameters</i>	65
3.3.5.	<i>Why to Use Pre-reaction Steps</i>	67

4.	Troubleshooting in Oven-based Systems	70
4.1.	Target Temperature Not Reached	70
4.2.	Significant Deviation of IR Temperatures	71
4.3.	Pressure Rate Increase Alarm	71
4.4.	Maintenance and Care of Sealed Vessels	73
4.5.	How to Check and Control Venting Activity in Closed Vessels with Venting Capability	74
4.5.1.	<i>How to Avoid Acid Venting</i>	75
4.5.2.	<i>Leakage</i>	75

Part B – Anton Paar Instrumentation and its Features

5.	Oven-based Systems	76
5.1.	Introduction to Oven-based Systems	76
5.1.1.	<i>Multiwave GO Plus</i>	77
5.1.2.	<i>Multiwave 5000</i>	78
5.2.	Rotor Overview	79
5.3.	Special Accessories for Multiwave 5000	83
5.3.1.	<i>Microwave Induced Combustion</i>	83
5.3.2.	<i>Disposable Glass Inserts</i>	84
5.3.3.	<i>Evaporation Accessory</i>	85
5.3.4.	<i>Drying Rotor</i>	86
5.4.	Applications in Multiwave GO Plus	87
5.4.1.	<i>Fat-rich Food Samples</i>	88
5.4.2.	<i>Organic Solvents & Alcoholic Beverages</i>	88
5.4.3.	<i>Polymer Samples</i>	89
5.4.4.	<i>Petroleum Samples</i>	89
5.4.5.	<i>Metals and Alloys</i>	90

6.	Autoclave Systems	91
6.1.	Microwave Digestion System Multiwave 7000	91
6.2.	Workflow – Differences from Oven-based Systems	92
6.3.	Racks and Vials	93
6.3.1.	<i>Pressure-sealed vials</i>	93
6.3.2.	<i>Sealed quartz vessels</i>	95
6.4.	Samples for Multiwave 7000	95
6.4.1.	<i>Samples with Unknown Reaction Behavior</i>	97
6.4.2.	<i>Increasing the Sample Weight</i>	97
6.4.3.	<i>Samples Containing Organic Solvents</i>	98
6.5.	Methods and Programs	98
6.5.1.	<i>How to Adapt Generic Methods</i>	99
7.	Appendix	101
7.1.	Abbreviations	101
7.2.	Glossary	102
7.3.	Application Matrix	103
7.4.	List of Common Standard Methods	105