

Semi-Solid and Solid Density Analyzers

Ultrapyc Series



Superior Density Measurements

**Achieve accurate and repeatable results
– no matter the sample.**

Gas pycnometry offers unmatched precision and repeatability for determining true and skeletal density – outperforming traditional techniques like liquid displacement or geometric methods. By using an inert gas to penetrate fine pores and irregular surfaces, gas pycnometers deliver highly reliable results, even in porous, irregular, or delicate solids.

From powders to slurries, this non-destructive method delivers consistent results across industries including pharmaceuticals, ceramics, polymers, and energy materials.

Taking density measurements to the next level

With superior features such as TruPyc dual reference chambers, the TruLock lid with repeatable closure, and an intuitive touchscreen interface, the Ultrapyc series makes density determination fast, accurate, and easy – bringing confidence to every measurement.



**Accommodates
samples from less
than 0.1 cm³ up to
135 cm³**

**Results in less
than a minute**

**Accurate densities
to within 0.015 %**

**Built-in and
customizable
methods**

**Peltier temperature
control
(3 °C to 60 °C)**

**Bi-directional
gas flow**



Pellets and monoliths

Reliably analyze a wide range of solid samples, from less than 0.1 cm³ to 135 cm³, to assess their internal porosity and structural integrity for quality assurance.



Powders

PowderProtect mode offers bi-directional gas expansion that eliminates the possibility of contaminating your instrument with fine powders.



Pastes and slurries

Peltier temperature control and disposable cups provide the flexibility to measure complex and hard-to-clean semi-solid samples with no mess.



Foams

Built-in methods compliant with ASTM D6226 determine a foam's open-cell content so you can predict and optimize their performance.

Four Solutions Tailored to Your Measurement Needs

| | Ultrapyc 3000 | Ultrapyc 5000 | Ultrapyc 7000 | Ultrapyc 7000 Micro |
|--|---------------|---------------|---------------|---------------------|
| TruPyc volume parity | ✓ | ✓ | ✓ | ✓ |
| TruLock repeatable sample chamber volume | ✓ | ✓ | ✓ | ✓ |
| PowderProtect bi-directional flow | | ✓ | ✓ | ✓ |
| Peltier temperature control | | ✓ | ✓ | ✓ |
| Method library | | | ✓ | ✓ |
| Foam mode | | | ✓ | |
| Small volume | | | | ✓ |



Ultrapyc 3000

Optimized volumes for accurate results
The Ultrapyc 3000 features TruPyc volume parity – two built-in reference chambers that automatically match your sample size for optimal accuracy. Its TruLock lid ensures repeatable chamber volume and best-in-class reproducibility across a wide range of sample volumes, from 0.1 cm³ to 135 cm³.

Ultrapyc 5000

Flexible control for complex samples
In addition to all Ultrapyc 3000 features, the Ultrapyc 5000 offers patented bi-directional control – choose reference-first (PowderProtect) or sample-first expansion to suit your material. With Peltier temperature control from 3 °C to 50 °C, the Ultrapyc 5000 ensures high-accuracy results even for volatile or temperature-sensitive samples.


Ultrapyc 7000

High throughput with advanced automation
The Ultrapyc 7000 builds on the Ultrapyc 5000 with an expanded temperature range up to 60 °C and a built-in method library for fast parameter switching – ideal for labs analyzing a variety of samples. It also includes dedicated foam methods (e.g., ASTM D6226) for streamlined testing of cellular materials.

Ultrapyc 7000 Micro

Precision for small samples
Designed for high-value or limited-quantity materials, the Ultrapyc 7000 Micro delivers all the benefits of the Ultrapyc 7000 while supporting ultra-small sample volumes – from less than 0.1 cm³, up to to 10 cm³ – without compromising accuracy or repeatability.

Streamline Your Experience



| SAMPLE NAME | DATE/TIME | INSTRUMENT TYPE | INSTRUMENT ALIAS | METHOD TYPE | STATUS |
|---------------------|-----------------------|-----------------|------------------|-------------|----------|
| 7/22/2025 (9 items) | | | | | |
| 1486 frac foam | 7/22/2025 12:05:33 PM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 frac foam | 7/22/2025 11:18:57 AM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 comp foam | 7/22/2025 10:44:51 AM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 comp foam | 7/22/2025 9:40:47 AM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 comp foam | 7/22/2025 9:18:13 AM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 frac foam | 7/22/2025 8:48:32 AM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 7/21/2025 (3 items) | | | | | |
| 1486 frac foam | 7/21/2025 6:12:39 PM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 comp foam | 7/21/2025 5:57:26 PM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 comp foam | 7/21/2025 5:36:59 PM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 unc foam | 7/21/2025 5:07:39 PM | Ultrapyc 7000 | TheKing | Measurement | ✓ |
| 1486 unc foam | 7/21/2025 10:52:30 AM | Ultrapyc 7000 | paul | Check | UNLOCKED |

1486 frac foam

Measurement 7/22/2025 11:18:57 AM 480003

Attachments ^

- 0_Ultrapyc_Report_20250722... PDF - 72 KB
- 0_Ultrapyc_Report_20250722... TXT - 2 KB

Measurement settings ^

Measurement setup ^

| | |
|------------------------|-----------------------|
| Sample ID | 1486 frac foam |
| Flow direction | Sample first |
| Thermal delay | Disabled |
| Start time | 7/22/2025 11:18:57 AM |
| Run time | 06:22:23.7342538 |
| Sample mass | 0.0010000 kg |
| Is balance sample mass | False |

Seamless data flow with AP Connect

Transfer your measurement data directly to Anton Paar's AP Connect lab execution system — no USBs, no paperwork, no delays. Results are securely stored, organized, and instantly available across your lab network, helping you streamline documentation, support compliance, and accelerate decision-making.

Optimize your workflow, eliminate errors

Connect your Ultrapyc to a balance and barcode reader for seamless data integration. Sample weights are transferred directly from the balance, while barcode scanning ensures fast, error-free sample identification. No manual entry means no transcription mistakes — just faster setup, cleaner records, and smoother lab operations.

Touchscreen simplicity, built-in intelligence

The Ultrapyc's 7-inch touchscreen and intuitive onboard software make operation effortless — from setting up methods to viewing results in real time. With guided workflows and multi-language support, users of all experience levels can run accurate measurements with confidence.



Find out more

Versatility Across Industries

The Ultrapyc series complies with many ASTM, ISO, and USP standard test methods used for determining the density of solids and semi-solids across a wide variety of industries.



Technical ceramics and refractories
Skeletal density reveals closed porosity and sintering efficiency, making it vital for structural ceramics and high-performance refractories. The Ultrapyc ensures reliable results for rigid and porous samples, helping you validate performance specifications and reduce material variability.



Food science and nutraceuticals
Understanding true density helps optimize texture, stability, and shelf life in powders, granulates, and compressed products. The Ultrapyc enables fast, non-destructive testing of moisture-sensitive or porous food materials, supporting product consistency and formulation improvements in R&D and QA settings.



Asphalt, bitumen, and road materials
Accurate skeletal density is essential for assessing compaction, void content, and durability in asphalt mixtures and fillers. The Ultrapyc delivers precise measurements even for viscous or semi-solid samples, while PowderProtect mode and Peltier temperature control ensure stability during testing.




Additive manufacturing and powder metallurgy
In additive manufacturing and powder metallurgy, knowing the true density of your powders ensures proper compaction, sintering behavior, and final part integrity. With Ultrapyc's PowderProtect mode and TruPyc volume parity, even fine, reactive powders are handled safely and accurately, without compromising data quality.




Pharmaceuticals
From active ingredients to tablet blends, accurate skeletal density is essential for controlling dosage uniformity, porosity, and dissolution rates. The Ultrapyc's small volume cells, intuitive touchscreen, and AP Connect integration streamline workflows and support data integrity.



Construction materials and binders
True density of cement, fly ash, and aggregates is critical for mix design, porosity evaluation, and long-term durability. The Ultrapyc accommodates large sample volumes up to 135 cm³ and offers robust, repeatable measurements, even in high-throughput environments.



Batteries and energy storage
For anodes, cathodes, and solid electrolytes, true density helps optimize energy density and performance. The Ultrapyc offers precise, temperature-controlled measurements and easy method switching for different chemistries – perfect for R&D teams characterizing multiple materials daily.



Polymers, plastics, and foams
Lightweight yet durable polymers and foams rely on accurate density profiling for quality assurance and material design. The Ultrapyc 7000 includes dedicated foam methods (e.g., ASTM D6226), while the touchscreen and guided workflows make testing quick and accessible – even for new operators.




Mining and oil exploration
True density is used to assess the purity and composition of the barite solids that control the density of drilling fluids, and to measure the volume of the solid phase of the concrete and cement materials used in well construction. Ultrapyc instruments can deliver true density results for powders, cores, and semi-solids in less than a minute.

| | Ultrapyc 3000 | Ultrapyc 5000 | Ultrapyc 7000 | Ultrapyc 7000 Micro |
|--|--|---------------------------------------|---|--|
| Performance specifications | | | | |
| Accuracy | 0.02 % ¹⁾ | | 0.015 % ¹⁾ | 0.075 % ²⁾ |
| Reference volumes | 2 (TruPyc) Nominally: 50 cm³, 8 cm³ | | | 2 (TruPyc) Nominally: 2 cm³, 5 cm³ |
| Inserts available | 135 cm³, 50 cm³, 10 cm³, 4.5 cm³, 1.8 cm³, 0.25 cm³ | | | 10 cm³, 4.5 cm³, 1.8 cm³, 0.25 cm³, 0.1 cm³ |
| Preparation modes | Flow, pulse, none | Flow, pulse, vacuum, none | | |
| Gas expansion direction | Sample chamber first | Bi-directional flow (PowderProtect) | | |
| Built-in temperature control | No | 3 °C to 50 °C ± 0.02 °C ³⁾ | 3 °C to 60 °C ± 0.02 °C ³⁾ | |
| Method library | No | | Yes | |
| Built-in foam methods and calculations | No | No | Yes | No |
| Sample chamber closure | Repeatable dual-turn self-aligning lid (TruLock) | | | |
| Display and control | 7-inch TFT WVGA (800 x 480 pixels); PCAP touchscreen | | | |
| Graphical user interface | Yes | | | |
| Transducer accuracy | Better than 0.1 % | | | |
| Connectivity | | | | |
| printer, bar/QR reader | Compatible via USB ports (4 total) | | | |
| PC / network connectivity | Yes (AP Connect) | | Yes (AP Connect or network shares) | |
| Language support | | | | |
| Languages | Chinese, English, French, German, Japanese, Korean Portuguese, Spanish, Turkish | | Chinese, English, French, German, Japanese, Korean, Polish, Portuguese, Spanish, Turkish | |
| Physical specifications | | | | |
| Width x depth x height | 27 cm (11 in) x 48 cm (19 in) x 25 cm (10 in) | | | |
| Weight | 10 kg (22 lbs) | | | |
| Gas used (not supplied) | Ultra-high purity helium, nitrogen, or any other non-reactive, non-corrosive gas (e.g., argon) regulated to no more than 30 psi/2.07 bar | | | |
| Power supply | External AC/DC adapter Input: AC 100 V to 240 V, 47 Hz to 63 Hz Output: 24 VDC, 3A | | External AC/DC adapter Input: 100 V to 240 V, 47 Hz to 63 Hz Output: 24 VDC, 5A | |
| Environmental specifications | | | | |
| Ambient temperature | 10 °C to 35 °C (50 °F to 95 °F) | | | |
| Air humidity | 10 % to 90 % RH non-condensing | | | |
| Altitude | Maximum 3,000 m (9,800 ft) | | | |
| Indoor / outdoor usage | Indoor only | | | |


1) With reference volume of 70.699 cm³
2) With reference volume of 2.145 cm³
3) For sample temperatures ≥ 15 °C, under standard laboratory conditions

Accessories


Foam sample preparation kit




Small-volume sample cells



Disposable cups



Non-elutriating cells



Selected international standards

| | | | |
|------------|--|-----------|------------------------------|
| ASTM B923 | Metal powders | ISO 787 | Pigments |
| ASTM C110 | Cements | ISO 4590 | Foam-rigid cellular plastics |
| ASTM C604 | Refractories | ISO 8130 | Coating powders - Part 2 |
| ASTM D2638 | Carbon (petroleum coke) | ISO 12154 | Solids |
| ASTM D4892 | Carbon (solid pitch) | USP 699 | Solids - pharmaceuticals |
| ASTM D5550 | Soil solids | | |
| ASTM D5965 | Coatings | | |
| ASTM D6093 | Pigments (clear or pigmented coatings) | | |
| ASTM D6226 | Rigids foams | | |



Our well-trained and certified technicians are ready to keep your instrument running smoothly.

Maximum uptime | Warranty program | Short response times | Global service network

