



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product Code(s)	AP-N7.5-100
Product Name	Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded
Form	Not applicable
Unique Formula Identifier (UFI)	A095-W03C-6002-8K08
Pure substance/mixture	Mixture

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Recommended use	Laboratory use
Uses advised against	No information available

### 1.3. Details of the supplier of the safety data sheet

#### Supplier

LGC Limited  
The Priestley Centre 10  
Priestley Road  
Surrey Research Park  
Guildford, England  
GU2 7XY  
Tel :+44 (0) 20 8943 7000  
Fax :+44 (0) 20 8943 2767  
eMail : gb@lgcstandards.com

Web : [www.lgcstandards.com](http://www.lgcstandards.com)

For further information, please contact

E-mail address sds-request@lgcstandards.com

### 1.4. Emergency telephone number

Emergency Telephone	For Hazardous Materials or Dangerous Goods Incident Spill, Leak, Fire Exposure, or Accident
---------------------	--



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

Call CHEMTREC:  
USA & Canada 1-800-424-9300  
Rest of the world +1 703-741-3877

Emergency Telephone - §45 - (EC)1272/2008	
Europe	112
Austria	No information available
Bulgaria	
Croatia	
Cyprus	
Czech Republic	
Denmark	
France	
Hungary	
Ireland	
Italy	
Lithuania	
Luxembourg	
Netherlands	
Norway	
Portugal	
Romania	
Slovakia	
Slovenia	
Spain	
Sweden	
Switzerland	

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to  
Regulation (EC) No. 1272/2008 [CLP]

Aspiration hazard	Category 1 - (H304)
-------------------	---------------------

### 2.2. Label elements

Contains 1-Decene, homopolymer, hydrogenated; Kerosine : distillates, hydrotreated light (D94/69 p1249) 4



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**



**Signal word**  
Danger

**Hazard statements**

H304 - May be fatal if swallowed and enters airways.

**Precautionary Statements - EU (§28, 1272/2008)**

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor

P331 - Do NOT induce vomiting

P405 - Store locked up

P501 - Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

**2.3. Other hazards**

**Other hazards** No information available.

**PBT or vPvB properties** The mixture does not contain any substances meeting the PBT or vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

**Endocrine Disruptor Information** This product does not contain any known or suspected endocrine disruptors.

Chemical name	EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances of Very High Concern (SVHC) for Authorisation	EU - REACH (1907/2006) - Endocrine Disruptor Assessment List of Substances
1-Decene, homopolymer, hydrogenated	-	-
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4	-	-

## SECTION 3: Composition/information on ingredients

**3.1 Substances**

Not applicable



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

## 3.2 Mixtures

**Chemical nature** Mixture of organic compounds.

Chemical name	Weight-%	REACH registration number	EC No. (Index No.)	Classification according to Regulation (EC) No. 1272/2008 [CLP]	Specific concentration limit (SCL)	M-Factor	M-Factor (long-term)	Notes
1-Decene, homopolymer, hydrogenated 68037-01-4	60 - 80	-	500-183-1	Asp. Tox. 1 (H304)	-	-	-	-
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8	40 - 60	-	265-149-8	Asp. Tox. 1 (H304)	-	-	-	-

*- Substances contained in this mixture with No registration number are under the REACH threshold in Article 6(1) and not subject to the registration requirements according to REACH Title II*

**Full text of H- and EUH-phrases: see section 16**

## Acute Toxicity Estimate

Chemical name	Oral LD50 mg/kg	Dermal LD50 mg/kg	Inhalation LC50 - 4 hour - dust/mist - mg/L	Inhalation LC50 - 4 hour - vapour - mg/L	Inhalation LC50 - 4 hour - gas - ppm
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8	5000	2000	No data available	No data available	No data available

This product does not contain candidate substances of very high concern at a concentration  $\geq 0.1\%$  (Regulation (EC) No. 1907/2006 (REACH), Article 59).

## SECTION 4: First aid measures

### 4.1. Description of first aid measures



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

---

<b>General advice</b>	Immediate medical attention is required. Show this safety data sheet to the doctor in attendance.
<b>Inhalation</b>	Aspiration into lungs can produce severe lung damage. If breathing has stopped, give artificial respiration. Get medical attention immediately. Remove to fresh air. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. If breathing is difficult, (trained personnel should) give oxygen. Delayed pulmonary edema may occur.
<b>Eye contact</b>	Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a doctor.
<b>Skin contact</b>	Wash skin with soap and water. In the case of skin irritation or allergic reactions see a doctor.
<b>Ingestion</b>	ASPIRATION HAZARD IF SWALLOWED - CAN ENTER LUNGS AND CAUSE DAMAGE. Do NOT induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration. Rinse mouth. Never give anything by mouth to an unconscious person. Get immediate medical attention.
<b>Self-protection of the first aider</b>	Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Use personal protective equipment as required.

#### **4.2. Most important symptoms and effects, both acute and delayed**

**Symptoms** Difficulty in breathing. Coughing and/ or wheezing. Dizziness.

**Effects of Exposure** No information available.

#### **4.3. Indication of any immediate medical attention and special treatment needed**

**Note to doctors** Because of the danger of aspiration, emesis or gastric lavage should not be used unless the risk is justified by the presence of additional toxic substances.

## **SECTION 5: Firefighting measures**

### **5.1. Extinguishing media**

**Suitable Extinguishing Media** Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

**Large Fire** CAUTION: Use of water spray when fighting fire may be inefficient.

**Unsuitable extinguishing media** Do not scatter spilled material with high pressure water streams.



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

---

## 5.2. Special hazards arising from the substance or mixture

**Specific hazards arising from the chemical** No information available.

## 5.3. Advice for firefighters

**Special protective equipment and precautions for fire-fighters** Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

## **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

**Personal precautions** Ensure adequate ventilation. Use personal protective equipment as required.

**Other information** Refer to protective measures listed in Sections 7 and 8.

**For emergency responders** Use personal protection recommended in Section 8.

### 6.2. Environmental precautions

**Environmental precautions** See Section 12 for additional Ecological Information.

### 6.3. Methods and material for containment and cleaning up

**Methods for containment** Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Dam up. Soak up with inert absorbent material. Pick up and transfer to properly labelled containers. Clean contaminated surface thoroughly.

**Prevention of secondary hazards** Clean contaminated objects and areas thoroughly observing environmental regulations.

### 6.4. Reference to other sections

**Reference to other sections** See section 8 for more information. See section 13 for more information.

## **SECTION 7: Handling and storage**

### 7.1. Precautions for safe handling

**Advice on safe handling** Ensure adequate ventilation.

---



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

**General hygiene considerations** Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Wash hands before breaks and after work. Wear suitable gloves and eye/face protection.

## 7.2. Conditions for safe storage, including any incompatibilities

**Storage Conditions** Please refer to the manufacturer's certificate for specific storage and transport temperature conditions. Store only in the original receptacle unless other advice is given on the CoA. Keep containers tightly closed in a dry, cool and well-ventilated place. Store locked up. Keep out of the reach of children. Store away from other materials.

**Storage class (TRGS 510)** LGK 10.

## 7.3. Specific end use(s)

**Risk Management Methods (RMM)** The information required is contained in this Safety Data Sheet.

## **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

#### **Exposure Limits**

Exposure Limits					
Chemical name	France	Germany TRGS	Germany DFG	Greece	Hungary
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8	-	TWA-AGW; ();	TWA-MAK: 5 mg/m³ ; II(4); respirable fraction TWA-MAK: 50 ppm; II(2); TWA-MAK: 350 mg/m³; II(2);	-	-
Chemical name	Sweden		Switzerland	United Kingdom	
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8	-		TWA-MAK: 50 ppm; vapour TWA-MAK: 350 mg/m³; vapour TWA-MAK: 5 mg/m³; aerosol, inhalable dust STEL-KZGW: 100 ppm; STEL-KZGW: 700 mg/m³;	-	

**Biological occupational exposure limits** This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

**Derived No Effect Level (DNEL) - Workers** No information available

**Derived No Effect Level (DNEL) - General Public** No information available.

**Predicted No Effect Concentration (PNEC)** No information available.

## 8.2. Exposure controls

### Personal protective equipment

**Eye/face protection** Avoid contact with eyes. Wear safety glasses with side shields (or goggles).

**Hand protection** The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374. Wear protective nitrile rubber gloves. Wear suitable gloves.

**Skin and body protection** Wear suitable protective clothing.

**Respiratory protection** Use appropriate respiratory protection. Consult with an industrial hygienist to determine the appropriate respiratory protection for your specific use of this material.

**Recommended filter type:** white. Particulates filter conforming to EN 143.

**Environmental exposure controls** Do not allow into any sewer, on the ground or into any body of water.

## **SECTION 9: Physical and chemical properties**

### 9.1. Information on basic physical and chemical properties

<b>Appearance</b>	Liquid
<b>Physical state</b>	Liquid
<b>Colour</b>	colourless
<b>Odour</b>	Characteristic
<b>Odour threshold</b>	No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
<b>Melting point / freezing point</b>	No data available	None known
<b>Boiling point or initial boiling point and boiling range</b>	No data available	None known





# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

Flammability	No data available	None known
Lower and upper explosion limit/flammability limit		None known
Lower explosion limit	No data available	
Upper explosion limit	No data available	
Flash point	No data available	None known
Autoignition temperature	No data available	None known
Decomposition temperature		None known
SADT (°C)	No data available	None known
pH	No data available	None known
pH (as aqueous solution)	No data available	None known
Kinematic viscosity	No data available	None known
Dynamic viscosity	No data available	None known
Solubility	No data available	None known
Water solubility	No data available	None known
Partition coefficient n-octanol/water (log value)	No data available	None known
Vapour pressure	No data available	None known
Density and/or relative density	No data available	None known
Bulk density	No data available	
Liquid Density	No data available	
Relative vapour density	No data available	None known
Particle characteristics		
Particle Size	No information available	
Particle Size Distribution	No information available	

## 9.2. Other information

### 9.2.1. Information with regards to physical hazard classes

No information available

### 9.2.2. Other safety characteristics

No information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

Reactivity No information available.

### 10.2. Chemical stability

Stability Stable under normal conditions.

Explosion data



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

---

**Sensitivity to mechanical impact** None.

**Sensitivity to static discharge** None.

## 10.3. Possibility of hazardous reactions

**Possibility of hazardous reactions** None under normal processing.

## 10.4. Conditions to avoid

**Conditions to avoid** None known based on information supplied.

## 10.5. Incompatible materials

**Incompatible materials** None known based on information supplied.

## 10.6. Hazardous decomposition products

**Hazardous decomposition products** None known based on information supplied.

## **SECTION 11: Toxicological information**

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Information on likely routes of exposure

##### **Product Information**

<b>Inhalation</b>	Specific test data for the substance or mixture is not available. Aspiration into lungs can produce severe lung damage. May cause pulmonary edema. Pulmonary edema can be fatal.
<b>Eye contact</b>	Specific test data for the substance or mixture is not available.
<b>Skin contact</b>	Specific test data for the substance or mixture is not available.
<b>Ingestion</b>	Specific test data for the substance or mixture is not available. Potential for aspiration if swallowed. May cause lung damage if swallowed. Aspiration may cause pulmonary edema and pneumonitis. May be fatal if swallowed and enters airways.

#### Symptoms related to the physical, chemical and toxicological characteristics

**Symptoms** Difficulty in breathing. Coughing and/ or wheezing. Dizziness.

**Acute toxicity** Based on available data, the classification criteria are not met.

---



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

## Numerical measures of toxicity

The following ATE values have been calculated for the mixture

ATEmix (oral)	99,999.00 mg/kg
ATEmix (dermal)	99,999.00 mg/kg
ATEmix (inhalation-gas)	99,999.00 ppm
ATEmix (inhalation-vapour)	99,999.00 mg/L
ATEmix (inhalation-dust/mist)	99,999.00 mg/L

## Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
1-Decene, homopolymer, hydrogenated	-	-	> 5.2 mg/L ( Rat ) 4 h
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4	> 5000 mg/kg ( Rat )	> 2000 mg/kg ( Rabbit )	> 5.2 mg/L ( Rat ) 4 h

## Delayed and immediate effects as well as chronic effects from short and long-term exposure

<b>Skin corrosion/irritation</b>	Based on available data, the classification criteria are not met.
<b>Serious eye damage/eye irritation</b>	Based on available data, the classification criteria are not met.
<b>Respiratory or skin sensitisation</b>	Based on available data, the classification criteria are not met.
<b>Germ cell mutagenicity</b>	Based on available data, the classification criteria are not met.
<b>Carcinogenicity</b>	Based on available data, the classification criteria are not met.
<b>Reproductive toxicity</b>	Based on available data, the classification criteria are not met.
<b>STOT - single exposure</b>	Based on available data, the classification criteria are not met.
<b>STOT - repeated exposure</b>	Based on available data, the classification criteria are not met.
<b>Aspiration hazard</b>	May be fatal if swallowed and enters airways.



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

## 11.2. Information on other hazards

### 11.2.1. Endocrine disrupting properties

**Endocrine disruption for human health** This product does not contain any known or suspected endocrine disruptors.

### 11.2.2. Other information

**Other adverse effects** No information available.

## **SECTION 12: Ecological information**

12.1. Toxicity Based on available data, the classification criteria are not met.

### **Aquatic toxicity**

#### **Component Information**

Chemical name	Fish	Crustacea	Algae/aquatic plants	Toxicity to microorganisms
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4	LC50: =2.2mg/L (96h, Lepomis macrochirus) LC50: =2.4mg/L (96h, Oncorhynchus mykiss) LC50: =45mg/L (96h, Pimephales promelas)	-	-	-

12.2. Persistence and degradability No information available.

### 12.3. Bioaccumulative potential

Chemical name	Partition coefficient	Bioconcentration factor (BCF)	Trophic magnification factor (TMF)
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4	-	159	-

12.4. Mobility in soil No information available.



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

## **12.5. Results of PBT and vPvB assessment**

This product does not contain any substances that are assessed to be a PBT or a vPvB.

Chemical name	PBT and vPvB assessment
1-Decene, homopolymer, hydrogenated	The substance is not PBT / vPvB
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4	The substance is not PBT / vPvB

## **12.6. Endocrine disrupting properties**

Based on available data, the classification criteria are not met.

## **12.7. Other adverse effects**

No information available.

## **PMT or vPvM properties**

Based on available data, the classification criteria are not met.

## **SECTION 13: Disposal considerations**

### **13.1. Waste treatment methods**

#### **Waste from residues/unused products**

Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

#### **Contaminated packaging**

Do not reuse empty containers.

## **SECTION 14: Transport information**

### **IATA**

14.1 UN number or ID number	Not regulated
14.2	
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable
14.6 Special precautions for user	
Special Provisions	None

### **IMDG**

14.1 UN number or ID number	Not regulated
14.2	
14.3 Transport hazard class(es)	Not regulated
14.4 Packing group	Not regulated
14.5 Environmental hazards	Not applicable



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

## 14.6 Special precautions for user

### Special Provisions

None

## 14.7 Maritime transport in bulk according to IMO instruments

No information available

## RID

### 14.1 UN number or ID number

Not regulated

### 14.2

### 14.3 Transport hazard class(es)

Not regulated

### 14.4 Packing group

Not regulated

### 14.5 Environmental hazards

Not applicable

## 14.6 Special precautions for user

### Special Provisions

None

## ADR

### 14.1 UN number or ID number

Not regulated

### 14.2 UN proper shipping name

Not regulated

### 14.3 Transport hazard class(es)

Not regulated

### 14.4 Packing group

Not regulated

### 14.5 Environmental hazards

Not applicable

## 14.6 Special precautions for user

### Special Provisions

None

## ADN

### 14.1 UN number or ID number

Not regulated

### 14.2

### 14.3 Transport hazard class(es)

Not regulated

### 14.4 Packing group

Not regulated

### 14.5 Environmental hazard

Not applicable

## 14.6 Special precautions for user

### Special Provisions

None

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### National regulations

#### France

#### Occupational Illnesses (R-463-3, France)

Chemical name	French RG number
Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8	RG 84

#### Germany



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

---

**Water hazard class (WGK)** slightly hazardous to water (WGK 1)

**Chemical Prohibition Ordinance (ChemVerbotsV)**  
Not applicable.

**TRGS 905** Not applicable

**Switzerland**

<b>Ordinance on the Incentive Tax on Volatile Organic Compounds (OVOC) SR 814.018</b>	Not applicable
<b>Storage of Hazardous Material</b>	Not applicable
<b>WPO (GSchV) SR 814.201; WPA (GSchG) SR 814.20</b>	Class B
<b>Major Accidents Ordinance SR 814.012</b>	Not applicable

**European Union**

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work.

**Authorisations and/or restrictions on use:**  
Use restricted. See item: 3.

**Persistent Organic Pollutants**  
Not applicable

**Ozone-depleting substances (ODS) regulation (EC) 2024/590**  
Not applicable.

**Explosives Precursors Marketing and Use (2019/1148)**  
Not applicable.

**International Inventories**

<b>TSCA</b>	U.S. INVENTORY (TSCA): Listed on inventory. For purposes of 40 CFR 720.36, this product is for Research and Development (R&D) Use Only
<b>DSL/NDSL</b>	Contact supplier for inventory compliance status
<b>EINECS/ELINCS</b>	Contact supplier for inventory compliance status

---



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

ENCS	Contact supplier for inventory compliance status
IECSC	Contact supplier for inventory compliance status
KECL	Contact supplier for inventory compliance status
PICCS	Contact supplier for inventory compliance status
AIIC	Contact supplier for inventory compliance status
NZIoC	Contact supplier for inventory compliance status
TCSI	Contact supplier for inventory compliance status

## Legend:

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory  
**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List  
**EINECS/ELINCS** - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances  
**ENCS** - Japan Existing and New Chemical Substances  
**IECSC** - China Inventory of Existing Chemical Substances  
**KECL** - Korean Existing Chemicals Inventory  
**PICCS** - Philippines Inventory of Chemicals and Chemical Substances  
**AICS** - Australian Inventory of Chemical Substances  
**NZIoC** - New Zealand Inventory of Chemicals  
**TCSI** - Taiwan Chemical Substance Inventory

## 15.2. Chemical safety assessment

**Chemical Safety Report** A Chemical Safety Assessment is not required for this substance

## SECTION 16: Other information

### Full text of any hazard and/or precautionary statements referred to under Sections 2-15

H304 - May be fatal if swallowed and enters airways  
P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor  
P331 - Do NOT induce vomiting  
P405 - Store locked up  
P501 - Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

### Key or legend to abbreviations and acronyms used in the safety data sheet

List may include phrases which are not applicable to this product

ACGIH	American Conference of Governmental Industrial Hygienists
AIDII	Italian Association of Industrial Hygienists
ADN	Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Europe)
ADR	Agreement concerning the International Carriage of Dangerous Goods by Road (Europe)
AIIC	Australian Inventory of Industrial Chemicals
ATE	Acute Toxicity Estimate





# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

ASTM	American Society for the Testing of Materials
bar	Biological Reference Values for Chemical Compounds in the Work Area
BAT	Biological tolerance values for occupational exposure
BEL	Biological exposure limits
bw	Body weight
Ceiling	Maximum limit value
CLP	Classification, Labelling and Packaging Regulation; Regulation (EC) No 1272/2008
CMR	Carcinogen, Mutagen or Reproductive Toxicant
DFG	German Research Foundation
DOT	Department of Transportation (United States)
DSL	Domestic Substances List (Canada)
ECHA	European Chemicals Agency
EC Number	European Community number
EmS	Emergency Schedule
ENCS	Existing and New Chemical Substances (Japan)
EPA	Environmental Protection Agency
EWC	European Waste Codes
GHS	Globally Harmonized System
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IBC	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
ICAO	International Civil Aviation Organisation
IECSC	Inventory of Existing Chemical Substances in China
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
ISO	International Organisation for Standardisation
KECI	Korean Existing Chemicals Inventory
LC50	Lethal Concentration to 50% of a test population
LD50	Lethal Dose to 50% of a test population (Median Lethal Dose)
MAK	Maximum Concentration at the Workplace
MAL	Measuring Technical Hygienic Air Needs
MARPOL	International Convention for the Prevention of Pollution from Ships
MDLPS	Ministry of Labour and Social Policy
n.o.s.	Not Otherwise Specified
NOAEC	No Observed Adverse Effect Concentration
NOAEL	No Observed Adverse Effect Level
NOELR	No Observable Effect Loading Rate
NZIoC	New Zealand Inventory of Chemicals
OECD	Organization for Economic Cooperation and Development
OEL	Occupational exposure limits
PBT	Persistent, Bioaccumulative and Toxic substance



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

PICCS	Philippines Inventory of Chemicals and Chemical Substances
PMT	Persistent, Mobile and Toxic
PPE	Personal protective equipment
QSAR	Quantitative Structure Activity Relationship
REACH	Registration, Evaluation, Authorisation, and Restriction of Chemicals (REACH) Regulation (EC 1907/2006)
RID	Agreement concerning the International Carriage of Dangerous Goods by Rail (Europe)
SADT	Self-Accelerating Decomposition Temperature
SAR	Structure-activity relationship
SDS	Safety Data Sheet
SL	Surface Limit
STEL	Short Term Exposure Limit
STOT RE	Specific target organ toxicity - Repeated exposure
STOT SE	Specific target organ toxicity - Single exposure
SVHC	Substance of very high concern
TCSI	Taiwan Chemical Substance Inventory
TDG	Transport of Dangerous Goods (Canada)
TRGS	Technical Rule for Hazardous Substances
TSCA	Toxic Substances Control Act (United States)
TWA	Time-Weighted Average
UN	United Nations
VOC	Volatile organic compounds
vPvB	Very Persistent and Very Bioaccumulative
vPvM	Very Persistent and Very Mobile
As	Allergenic substance
C	Carcinogen
DS	Dermal Sensitizer
Ot	Ototoxicant
pOt	Ototoxicant - potential to cause hearing disorders
PS	Photosensitizer
RS	Respiratory Sensitizer
S	Sensitizer
poS	Sensitizer - capable of causing occupational asthma
Sa	Simple asphyxiant
Sd	Skin designation
pSd	Skin designation - potential for cutaneous absorption
Sdv	Skin designation - vacated
Sk	Skin notation
dSk	Skin notation - danger of cutaneous absorption
pSk	Skin notation - potential for cutaneous absorption

Classification procedure



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

Classification according to Regulation (EC) No. 1272/2008 [CLP]	Method Used
Acute oral toxicity	Calculation method
Acute dermal toxicity	Calculation method
Acute inhalation toxicity - gas	Calculation method
Acute inhalation toxicity - vapour	Calculation method
Acute inhalation toxicity - dust/mist	Calculation method
Skin corrosion/irritation	Calculation method
Serious eye damage/eye irritation	Calculation method
Respiratory sensitisation	Calculation method
Skin sensitisation	Calculation method
Mutagenicity	Calculation method
Carcinogenicity	Calculation method
Reproductive toxicity	Calculation method
STOT - single exposure	Calculation method
STOT - repeated exposure	Calculation method
Chronic aquatic toxicity	Calculation method
Acute aquatic toxicity	Calculation method
Aspiration hazard	Calculation method
Ozone	Calculation method

## Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)  
U.S. Environmental Protection Agency ChemView Database  
European Food Safety Authority (EFSA)  
European Chemicals Agency (ECHA) Committee for Risk Assessment (ECHA\_RAC)  
European Chemicals Agency (ECHA) (ECHA\_API)  
Environmental Protection Agency  
Acute Exposure Guideline Level(s) (AEGL(s))  
U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act  
U.S. Environmental Protection Agency High Production Volume Chemicals  
Food Research Journal  
Hazardous Substance Database  
International Uniform Chemical Information Database (IUCLID)  
National Institute of Technology and Evaluation (NITE)  
Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)  
NIOSH (National Institute for Occupational Safety and Health)  
National Library of Medicine's ChemID Plus (NLM CIP)  
National Library of Medicine's PubMed database (NLM PUBMED)  
U.S. National Toxicology Program (NTP)  
New Zealand's Chemical Classification and Information Database (CCID)  
Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications  
Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme  
Organisation for Economic Co-operation and Development Screening Information Data Set  
World Health Organization



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 27-Jan-2026

Revision Number 1.02

**AP-N7.5-100 - Viscosity Reference Standard General Purpose Type N7.5, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded**

---

Revision date 27-Jan-2026

## Disclaimer

The information in this safety data sheet (SDS) has been prepared with due care and is true and accurate to the best of our knowledge. The user must determine the suitability of the information for its particular purpose, ensure compliance with existing laws and regulations, and be aware that other or additional safety or performance considerations may arise when using, handling and/ or storing the material. The information in this SDS does not purport to be all inclusive or a guarantee as to the properties of the material supplied, and should be used only as a guide. LGC makes no warranties or representations as to the accuracy and completeness of the information contained herein, shall not be held responsible for the suitability of this information for the user's intended purposes or the consequences of such use, and shall not be liable for any damage or loss, howsoever arising, direct or otherwise.

**End of Safety Data Sheet**