



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

| | |
|---------------------------------|---|
| Product Code(s) | AP-S3-100 |
| Product Name | Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded |
| Form | Not applicable |
| Unique Formula Identifier (UFI) | HJ25-X0AK-600X-NY9E |
| Pure substance/mixture | Mixture |

1.2. Relevant identified uses of the substance or mixture and uses advised against

| | |
|----------------------|--------------------------|
| Recommended use | Laboratory use |
| Uses advised against | No information available |

1.3. Details of the supplier of the safety data sheet

Supplier

LGC Limited
Queens Road
Teddington
Middlesex TW11 0LY
UNITED KINGDOM
:+44 (0) 20 8943 7000
Fax :+44 (0) 20 8943 2767
eMail : gb@lgcstandards.com

Web : www.lgcstandards.com

For further information, please contact

E-mail address sds-request@lgcstandards.com

1.4. Emergency telephone number

Emergency Telephone For Hazardous Materials or Dangerous Goods Incident
Spill, Leak, Fire Exposure, or Accident
Call CHEMTREC:



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

USA & Canada 1-800-424-9300
Rest of the world +1 703-741-3877

| Emergency Telephone - §45 - (EC)1272/2008 | |
|---|--------------------------|
| Europe | 112 |
| Austria | No information available |
| Bulgaria | |
| Croatia | |
| Cyprus | |
| Czech Republic | |
| Denmark | |
| France | |
| Hungary | |
| Ireland | |
| Italy | |
| Lithuania | |
| Luxembourg | |
| Netherlands | |
| Norway | |
| Portugal | |
| Romania | |
| Slovakia | |
| Slovenia | |
| Spain | |
| Sweden | |
| Switzerland | |

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to
Regulation (EC) No. 1272/2008 [CLP]

| | |
|-------------------|---------------------|
| Aspiration hazard | Category 1 - (H304) |
|-------------------|---------------------|

2.2. Label elements

Contains Kerosine : distillates, hydrotreated light (D94/69 p1249) 4; White mineral oil



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded



Signal word
Danger

Hazard statements

H304 - May be fatal if swallowed and enters airways.

Precautionary Statements - EU (§28, 1272/2008)

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor

P331 - Do NOT induce vomiting

P405 - Store locked up

P501 - Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

2.3. Other hazards

No information available.

PBT & vPvB

This mixture contains no substance considered to be persistent, bioaccumulating or toxic (PBT). This mixture contains no substance considered to be very persistent nor very bioaccumulating (vPvB).

Endocrine Disruptor Information This product does not contain any known or suspected endocrine disruptors.

| Chemical name | EU - REACH (1907/2006) - Article 59(1) - Candidate List of Substances of Very High Concern (SVHC) for Authorisation | EU - REACH (1907/2006) - Endocrine Disruptor Assessment List of Substances |
|---|---|--|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 | - | - |
| White mineral oil | - | - |

SECTION 3: Composition/information on ingredients

3.1 Substances

Not applicable



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

3.2 Mixtures

Chemical nature Mixture of organic compounds.

| Chemical name | Weight-% | REACH registration number | EC No (EU Index No) | Classification according to Regulation (EC) No. 1272/2008 [CLP] | Specific concentration limit (SCL) | M-Factor | M-Factor (long-term) | Notes |
|--|----------|---------------------------|---------------------|---|------------------------------------|----------|----------------------|-------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8 | 80 - 100 | - | 265-149-8 | Asp. Tox. 1 (H304) | - | - | - | - |
| White mineral oil 8042-47-5 | 20 - 40 | - | 232-455-8 | Asp. Tox. 1 (H304) | - | - | - | - |

- Substances contained in this mixture with No registration number are under the REACH threshold in Article 6(1) and not subject to the registration requirements according to REACH Title II

Full text of H- and EUH-phrases: see section 16

Acute Toxicity Estimate

If LD50/LC50 data is not available or does not correspond to the classification category, then the appropriate conversion value from CLP Annex I, Table 3.1.2, is used to calculate the acute toxicity estimate (ATEmix) for classifying a mixture based on its components

| Chemical name | Oral LD50 mg/kg | Dermal LD50 mg/kg | Inhalation LC50 - 4 hour - dust/mist - mg/L | Inhalation LC50 - 4 hour - vapour - mg/L | Inhalation LC50 - 4 hour - gas - ppm |
|--|-----------------|-------------------|---|--|--------------------------------------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8 | 5000 | 2000 | No data available | No data available | No data available |
| White mineral oil 8042-47-5 | 5000 | No data available | No data available | No data available | No data available |

This product does not contain candidate substances of very high concern at a concentration $\geq 0.1\%$ (Regulation (EC) No. 1907/2006 (REACH), Article 59).

SECTION 4: First aid measures

4.1. Description of first aid measures



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | |
|---|--|
| General advice | Immediate medical attention is required. Show this safety data sheet to the doctor in attendance. |
| Inhalation | Aspiration into lungs can produce severe lung damage. If breathing has stopped, give artificial respiration. Get medical attention immediately. Remove to fresh air. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. If breathing is difficult, (trained personnel should) give oxygen. Delayed pulmonary edema may occur. |
| Eye contact | Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a doctor. |
| Skin contact | Wash skin with soap and water. In the case of skin irritation or allergic reactions see a doctor. |
| Ingestion | ASPIRATION HAZARD IF SWALLOWED - CAN ENTER LUNGS AND CAUSE DAMAGE. Do NOT induce vomiting. If vomiting occurs spontaneously, keep head below hips to prevent aspiration. Rinse mouth. Never give anything by mouth to an unconscious person. Get immediate medical attention. |
| Self-protection of the first aider | Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Use personal protective equipment as required. |

4.2. Most important symptoms and effects, both acute and delayed

Symptoms Difficulty in breathing. Coughing and/ or wheezing. Dizziness.

Effects of Exposure No information available.

4.3. Indication of any immediate medical attention and special treatment needed

Note to doctors Because of the danger of aspiration, emesis or gastric lavage should not be used unless the risk is justified by the presence of additional toxic substances.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable Extinguishing Media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Large Fire CAUTION: Use of water spray when fighting fire may be inefficient.

Unsuitable extinguishing media Do not scatter spilled material with high pressure water streams.



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

5.2. Special hazards arising from the substance or mixture

Specific hazards arising from the chemical No information available.

5.3. Advice for firefighters

Special protective equipment and precautions for fire-fighters Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation. Use personal protective equipment as required.

Other information Refer to protective measures listed in Sections 7 and 8.

For emergency responders Use personal protection recommended in Section 8.

6.2. Environmental precautions

Environmental precautions See Section 12 for additional Ecological Information.

6.3. Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Take up mechanically, placing in appropriate containers for disposal.

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

6.4. Reference to other sections

Reference to other sections See section 8 for more information. See section 13 for more information.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling Ensure adequate ventilation.



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

General hygiene considerations Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Wash hands before breaks and after work. Wear suitable gloves and eye/face protection.

7.2. Conditions for safe storage, including any incompatibilities

Storage Conditions Please refer to the manufacturer's certificate for specific storage and transport temperature conditions. Store only in the original receptacle unless other advice is given on the CoA. Keep containers tightly closed in a dry, cool and well-ventilated place. Store locked up. Keep out of the reach of children. Store away from other materials.

Storage class (TRGS 510) LGK 10.

7.3. Specific end use(s)

Risk Management Methods (RMM) The information required is contained in this Safety Data Sheet.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Exposure Limits

| Chemical name | France | Germany TRGS | Germany DFG | Greece | Hungary |
|---|----------|--------------------------|---|--|--------------------------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8 | - | - | TWA: 5 mg/m ³ TWA: 50 ppm TWA: 350 mg/m ³ Peak: 20 mg/m ³ Peak: 100 ppm Peak: 700 mg/m ³ | - | - |
| White mineral oil 8042-47-5 | - | TWA: 5 mg/m ³ | TWA: 5 mg/m ³ Peak: 20 mg/m ³ | - | TWA: 5 mg/m ³ |
| Chemical name | Ireland | Italy MDLPS | Italy AIDII | Latvia | Lithuania |
| White mineral oil 8042-47-5 | - | - | - | TWA: 5 mg/m ³ | - |
| Chemical name | Portugal | Romania | Slovakia | Slovenia | Spain |
| White mineral oil 8042-47-5 | - | - | - | TWA: 5 mg/m ³ STEL: 20 mg/m ³ | - |
| Chemical name | Sweden | | Switzerland | | United Kingdom |
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8 | - | | TWA: 50 ppm TWA: 350 mg/m ³ TWA: 5 mg/m ³ | | - |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | | | |
|--------------------------------|---|--|---|
| | | STEL: 100 ppm STEL: 700 mg/m ³ | |
| White mineral oil 8042-47-5 | - | TWA: 5 mg/m ³ | - |

Biological occupational exposure limits This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

Derived No Effect Level (DNEL) - Workers

| Chemical name | Oral | Dermal | Inhalation |
|--------------------------------|------|-----------------------------|----------------------------------|
| White mineral oil 8042-47-5 | - | 217.05 mg/kg bw/day [4] [6] | 164.56 mg/m ³ [4] [6] |

Notes

[4] Systemic health effects.
[6] Long term.

Derived No Effect Level (DNEL) - General Public

| Chemical name | Oral | Dermal | Inhalation |
|--------------------------------|-------------------------|--------|---------------------------------|
| White mineral oil 8042-47-5 | 25 mg/kg bw/day [4] [6] | - | 34.78 mg/m ³ [4] [6] |

Notes

[4] Systemic health effects.
[6] Long term.

Predicted No Effect Concentration (PNEC) No information available.

8.2. Exposure controls

Personal protective equipment

Eye/face protection Avoid contact with eyes. Wear safety glasses with side shields (or goggles). If splashes are likely to occur, wear safety glasses with side-shields.

Hand protection The protective gloves to be used must comply with the specifications of EC Directive 89/686/EEC and the related standard EN374. Wear protective nitrile rubber gloves.



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | |
|--|---|
| Skin and body protection | Wear suitable protective clothing. |
| Respiratory protection | No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required. |
| General hygiene considerations | Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Wash hands before breaks and after work. Wear suitable gloves and eye/face protection. |
| Environmental exposure controls | Do not allow into any sewer, on the ground or into any body of water. |

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| | |
|------------------------|--------------------------|
| Appearance | Liquid |
| Physical state | Liquid |
| Colour | colourless |
| Odour | Characteristic |
| Odour threshold | No information available |

| Property | Values | Remarks • Method |
|--|-------------------|-------------------------|
| Melting point / freezing point | No data available | None known |
| Initial boiling point and boiling range | No data available | None known |
| Flammability | No data available | None known |
| Flammability Limit in Air | | None known |
| Upper flammability or explosive limits | No data available | |
| Lower flammability or explosive limits | No data available | |
| Flash point | No data available | None known |
| Autoignition temperature | No data available | None known |
| Decomposition temperature | | None known |
| SADT (°C) | No data available | None known |
| pH | No data available | None known |
| pH (as aqueous solution) | No data available | None known |
| Kinematic viscosity | No data available | None known |
| Dynamic viscosity | No data available | None known |
| Water solubility | No data available | None known |
| Solubility(ies) | No data available | None known |
| Partition coefficient | No data available | None known |
| Vapour pressure | No data available | None known |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | | |
|----------------------------|--------------------------|------------|
| Relative density | No data available | None known |
| Bulk density | No data available | |
| Liquid Density | No data available | |
| Relative vapour density | No data available | None known |
| Particle characteristics | | |
| Particle Size | No information available | |
| Particle Size Distribution | No information available | |

9.2. Other information

9.2.1. Information with regards to physical hazard classes

No information available

9.2.2. Other safety characteristics

No information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Reactivity No information available.

10.2. Chemical stability

Stability Stable under normal conditions.

Explosion data

Sensitivity to mechanical impact None.

Sensitivity to static discharge None.

10.3. Possibility of hazardous reactions

Possibility of hazardous reactions None under normal processing.

10.4. Conditions to avoid

Conditions to avoid None known based on information supplied.

10.5. Incompatible materials

Incompatible materials None known based on information supplied.

10.6. Hazardous decomposition products

Hazardous decomposition products None known based on information supplied.



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Information on likely routes of exposure

| | |
|--------------|--|
| Inhalation | Specific test data for the substance or mixture is not available. Aspiration into lungs can produce severe lung damage. May cause pulmonary edema. Pulmonary edema can be fatal. |
| Eye contact | Specific test data for the substance or mixture is not available. |
| Skin contact | Specific test data for the substance or mixture is not available. |
| Ingestion | Specific test data for the substance or mixture is not available. Potential for aspiration if swallowed. May cause lung damage if swallowed. Aspiration may cause pulmonary edema and pneumonitis. May be fatal if swallowed and enters airways. |

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms Difficulty in breathing. Coughing and/ or wheezing. Dizziness.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Acute toxicity Based on available data, the classification criteria are not met.

Numerical measures of toxicity

The following values are calculated based on chapter 3.1 of the GHS document:

ATEmix (oral) 99,999.00 mg/kg
ATEmix (dermal) 99,999.00 mg/kg
ATEmix (inhalation-gas) 99,999.00 ppm
ATEmix (inhalation-vapour) 99,999.00 mg/l
ATEmix (inhalation-dust/mist) 99,999.00 mg/l

| Chemical name | Oral LD50 | Dermal LD50 | Inhalation LC50 |
|---|----------------------|-------------------------|------------------------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 | > 5000 mg/kg (Rat) | > 2000 mg/kg (Rabbit) | > 5.2 mg/L (Rat) 4 h |
| White mineral oil | > 5000 mg/kg (Rat) | > 2000 mg/kg (Rabbit) | - |

Skin corrosion/irritation Based on available data, the classification criteria are not met.

Serious eye damage/eye irritation Based on available data, the classification criteria are not met.

Respiratory or skin sensitisation Based on available data, the classification criteria are not met.



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

Germ cell mutagenicity Based on available data, the classification criteria are not met.

Carcinogenicity Based on available data, the classification criteria are not met.

Reproductive toxicity Based on available data, the classification criteria are not met.

STOT - single exposure Based on available data, the classification criteria are not met.

STOT - repeated exposure Based on available data, the classification criteria are not met.

Aspiration hazard May be fatal if swallowed and enters airways.

11.2. Information on other hazards

11.2.1. Endocrine disrupting properties

Endocrine disrupting properties This product does not contain any known or suspected endocrine disruptors.

11.2.2. Other information

Other adverse effects No information available.

SECTION 12: Ecological information

12.1. Toxicity

Ecotoxicity

The environmental impact of this product has not been fully investigated.

| Chemical name | Algae/aquatic plants | Fish | Toxicity to microorganisms | Crustacea |
|---|----------------------|---|----------------------------|-----------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 | - | LC50: =2.2mg/L (96h, Lepomis macrochirus) LC50: =2.4mg/L (96h, Oncorhynchus mykiss) LC50: =45mg/L (96h, | - | - |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | | | | |
|-------------------|---|---|---|---|
| | | Pimephales promelas | | |
| White mineral oil | - | LC50: >10000mg/L (96h, Lepomis macrochirus) | - | - |

12.2. Persistence and degradability

Persistence and degradability No information available.

12.3. Bioaccumulative potential

Bioaccumulation

| Chemical name | Partition coefficient |
|-------------------|-----------------------|
| White mineral oil | 6 |

12.4. Mobility in soil

Mobility in soil No information available.

12.5. Results of PBT and vPvB assessment

PBT and vPvB assessment Based on available data, the classification criteria are not met.

| Chemical name | PBT and vPvB assessment |
|---|---------------------------------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 | The substance is not PBT / vPvB |
| White mineral oil | The substance is not PBT / vPvB |

12.6. Endocrine disrupting properties

Endocrine disrupting properties Based on available data, the classification criteria are not met.

12.7. Other adverse effects

Other adverse effects No information available.

PMT or vPvM properties Based on available data, the classification criteria are not met.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste from residues/unused products Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

Contaminated packaging Do not reuse empty containers.

SECTION 14: Transport information



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

IATA

| | |
|-----------------------------------|----------------|
| 14.1 UN number or ID number | Not regulated |
| 14.2 UN proper shipping name | Not regulated |
| 14.3 Transport hazard class(es) | Not regulated |
| 14.4 Packing group | Not regulated |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | None |

IMDG

| | |
|--|--------------------------|
| 14.1 UN number or ID number | Not regulated |
| 14.2 UN proper shipping name | Not regulated |
| 14.3 Transport hazard class(es) | Not regulated |
| 14.4 Packing group | Not regulated |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | None |
| 14.7 Maritime transport in bulk according to IMO instruments | No information available |

RID

| | |
|-----------------------------------|----------------|
| 14.1 UN number or ID number | Not regulated |
| 14.2 UN proper shipping name | Not regulated |
| 14.3 Transport hazard class(es) | Not regulated |
| 14.4 Packing group | Not regulated |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | None |

ADR

| | |
|-----------------------------------|----------------|
| 14.1 UN number or ID number | Not regulated |
| 14.2 UN proper shipping name | Not regulated |
| 14.3 Transport hazard class(es) | Not regulated |
| 14.4 Packing group | Not regulated |
| 14.5 Environmental hazards | Not applicable |
| 14.6 Special precautions for user | |
| Special Provisions | None |

ADN

| | |
|---------------------------------|----------------|
| 14.1 UN number or ID number | Not regulated |
| 14.2 UN proper shipping name | Not regulated |
| 14.3 Transport hazard class(es) | Not regulated |
| 14.4 Packing group | Not regulated |
| 14.5 Environmental hazard | Not applicable |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

14.6 Special precautions for user
Special Provisions None

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulations

France

Occupational Illnesses (R-463-3, France)

| Chemical name | French RG number | Title |
|---|------------------|-------|
| Kerosine : distillates, hydrotreated light (D94/69 p1249) 4 64742-47-8 | RG 84 | - |
| White mineral oil 8042-47-5 | RG 36bis | - |

Germany

Water hazard class (WGK) slightly hazardous to water (WGK 1)
Chemical Prohibition Ordinance (ChemVerbotsV) Not applicable

TRGS 905 Not applicable

Switzerland

Ordinance on the Incentive Tax on Volatile Organic Compounds (OVOC) SR 814.018 Not applicable
Storage of Hazardous Material Not applicable
WPO (GSchV) SR 814.201; WPA (GSchG) SR 814.20 Class B
Major Accidents Ordinance SR 814.012 Not applicable

European Union

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work.

Authorisations and/or restrictions on use:

This product does not contain substances subject to authorisation (Regulation (EC) No. 1907/2006 (REACH), Annex XIV) This product does not contain substances subject to restriction (Regulation (EC) No. 1907/2006 (REACH), Annex XVII)

DIRECTIVE (EU) 2021/1187 on the marketing and use of explosives precursors

Not applicable



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

Persistent Organic Pollutants
Not applicable

Ozone-depleting substances (ODS) regulation (EC) 2024/590
Not applicable

EU - Plant Protection Products (1107/2009/EC)

| Chemical name | EU - Plant Protection Products (1107/2009/EC) |
|-------------------------------|---|
| White mineral oil - 8042-47-5 | Plant protection agent |

International Inventories

TSCA

LGC, to the best of its ability, has confirmed that the chemical substances in this product are listed as "Active" in the EPA (Environmental Protection Agency) "TSCA Inventory Notification (Active-Inactive) Requirements Rule" ("the Final Rule") of Feb 2019, as amended Feb 2021."

DSL/NDSL

Contact supplier for inventory compliance status

EINECS/ELINCS

Contact supplier for inventory compliance status

ENCS

Contact supplier for inventory compliance status

IECSC

Contact supplier for inventory compliance status

KECL

Contact supplier for inventory compliance status

PICCS

Contact supplier for inventory compliance status

AIIC

Contact supplier for inventory compliance status

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing Chemicals Inventory

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AIIC - Australian Inventory of Chemical Substances

15.2. Chemical safety assessment

Chemical Safety Report

A Chemical Safety Assessment is not required for this substance



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

SECTION 16: Other information

Key or legend to abbreviations and acronyms used in the safety data sheet

Full text of any hazard and/or precautionary statements referred to under Sections 2-15

H304 - May be fatal if swallowed and enters airways

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor

P331 - Do NOT induce vomiting

P405 - Store locked up

P501 - Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

Legend

| | |
|-----------|---|
| ACGIH | American Conference of Governmental Industrial Hygienists |
| AIDII | Italian Association of Industrial Hygienists |
| ADN | Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Europe) |
| ADR | Agreement concerning the International Carriage of Dangerous Goods by Road (Europe) |
| AIIC | Australian Inventory of Industrial Chemicals |
| ATE | Acute Toxicity Estimate |
| ASTM | American Society for the Testing of Materials |
| bar | Biological Reference Values for Chemical Compounds in the Work Area |
| BAT | Biological tolerance values for occupational exposure |
| BEL | Biological exposure limits |
| bw | Body weight |
| Ceiling | Maximum limit value |
| CLP | Classification, Labelling and Packaging Regulation; Regulation (EC) No 1272/2008 |
| CMR | Carcinogen, Mutagen or Reproductive Toxicant |
| DFG | German Research Foundation |
| DOT | Department of Transportation (United States) |
| DSL | Domestic Substances List (Canada) |
| ECHA | European Chemicals Agency |
| EC Number | European Community number |
| EmS | Emergency Schedule |
| ENCS | Existing and New Chemical Substances (Japan) |
| EPA | Environmental Protection Agency |
| EWC | European Waste Codes |
| GHS | Globally Harmonized System |
| IARC | International Agency for Research on Cancer |
| IATA | International Air Transport Association |
| IBC | International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk |
| ICAO | International Civil Aviation Organisation |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | |
|---------|---|
| IECSC | Inventory of Existing Chemical Substances in China |
| IMDG | International Maritime Dangerous Goods |
| IMO | International Maritime Organization |
| ISO | International Organisation for Standardisation |
| KECI | Korean Existing Chemicals Inventory |
| LC50 | Lethal Concentration to 50% of a test population |
| LD50 | Lethal Dose to 50% of a test population (Median Lethal Dose) |
| MAL | Measuring Technical Hygienic Air Needs |
| MARPOL | International Convention for the Prevention of Pollution from Ships |
| MDLPS | Ministry of Labour and Social Policy |
| n.o.s. | Not Otherwise Specified |
| NOAEC | No Observed Adverse Effect Concentration |
| NOAEL | No Observed Adverse Effect Level |
| NOELR | No Observable Effect Loading Rate |
| NZIoC | New Zealand Inventory of Chemicals |
| OECD | Organization for Economic Cooperation and Development |
| OEL | Occupational exposure limits |
| PBT | Persistent, Bioaccumulative and Toxic substance |
| PICCS | Philippines Inventory of Chemicals and Chemical Substances |
| PMT | Persistent, Mobile and Toxic |
| PPE | Personal protective equipment |
| QSAR | Quantitative Structure Activity Relationship |
| REACH | Registration, Evaluation, Authorisation, and Restriction of Chemicals (REACH) Regulation (EC 1907/2006) |
| RID | Agreement concerning the International Carriage of Dangerous Goods by Rail (Europe) |
| SADT | Self-Accelerating Decomposition Temperature |
| SAR | Structure-activity relationship |
| SDS | Safety Data Sheet |
| SL | Surface Limit |
| STEL | Short Term Exposure Limit |
| STOT RE | Specific target organ toxicity - Repeated exposure |
| STOT SE | Specific target organ toxicity - Single exposure |
| SVHC | Substance of very high concern |
| TCSI | Taiwan Chemical Substance Inventory |
| TDG | Transport of Dangerous Goods (Canada) |
| TRGS | Technical Rule for Hazardous Substances |
| TSCA | Toxic Substances Control Act (United States) |
| TWA | Time-Weighted Average |
| UN | United Nations |
| VOC | Volatile organic compounds |
| vPvB | Very Persistent and Very Bioaccumulative |
| vPvM | Very Persistent and Very Mobile |



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

| | |
|------|--------------------|
| Sen+ | Sensitiser |
| Sk* | Skin designation |
| ** | Hazard Designation |

| Classification procedure | |
|---|--------------------|
| Classification according to Regulation (EC) No. 1272/2008 [CLP] | Method Used |
| Acute oral toxicity | Calculation method |
| Acute dermal toxicity | Calculation method |
| Acute inhalation toxicity - gas | Calculation method |
| Acute inhalation toxicity - vapour | Calculation method |
| Acute inhalation toxicity - dust/mist | Calculation method |
| Skin corrosion/irritation | Calculation method |
| Serious eye damage/eye irritation | Calculation method |
| Respiratory sensitisation | Calculation method |
| Skin sensitisation | Calculation method |
| Mutagenicity | Calculation method |
| Carcinogenicity | Calculation method |
| Reproductive toxicity | Calculation method |
| STOT - single exposure | Calculation method |
| STOT - repeated exposure | Calculation method |
| Chronic aquatic toxicity | Calculation method |
| Acute aquatic toxicity | Calculation method |
| Aspiration hazard | Calculation method |
| Ozone | Calculation method |

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)
U.S. Environmental Protection Agency ChemView Database
European Food Safety Authority (EFSA)
European Chemicals Agency (ECHA) Committee for Risk Assessment (ECHA_RAC)
European Chemicals Agency (ECHA) (ECHA_API)
Environmental Protection Agency
Acute Exposure Guideline Level(s) (AEGL(s))
U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
U.S. Environmental Protection Agency High Production Volume Chemicals
Food Research Journal
Hazardous Substance Database
International Uniform Chemical Information Database (IUCLID)
National Institute of Technology and Evaluation (NITE)
Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)
NIOSH (National Institute for Occupational Safety and Health)
National Library of Medicine's ChemID Plus (NLM CIP)
National Library of Medicine's PubMed database (NLM PUBMED)



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Regulation (EC) No. 1907/2006 and Regulation (EC) No. 1272/2008

Revision date 09-Jun-2025

Revision Number 1.02

AP-S3-100 - Viscosity Reference Standard General Purpose Type S3, UKAS ISO 17025 / ISO 17034, 100 mL / AP Dual-branded

U.S. National Toxicology Program (NTP)

New Zealand's Chemical Classification and Information Database (CCID)

Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications

Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme

Organisation for Economic Co-operation and Development Screening Information Data Set

World Health Organization

Revision date 09-Jun-2025

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006

Disclaimer

The information in this safety data sheet (SDS) has been prepared with due care and is true and accurate to the best of our knowledge. The user must determine the suitability of the information for its particular purpose, ensure compliance with existing laws and regulations, and be aware that other or additional safety or performance considerations may arise when using, handling and/ or storing the material. The information in this SDS does not purport to be all inclusive or a guarantee as to the properties of the material supplied, and should be used only as a guide. LGC makes no warranties or representations as to the accuracy and completeness of the information contained herein, shall not be held responsible for the suitability of this information for the user's intended purposes or the consequences of such use, and shall not be liable for any damage or loss, howsoever arising, direct or otherwise.

End of Safety Data Sheet