

Please find enclosed a list of already tested printers:

USB connection	Ethernet connection
HP Color Laserjet CP4025	HP Business Inkjet 2300
HP Deskjet 5560	HP Color Laserjet 3800n
HP Deskjet 6940	HP Color Laserjet M553
HP Deskjet 6980	HP Laserjet 4050N
HP Deskjet 6983	HP Laserjet 4100N
HP Color Laserjet M553	HP Laserjet 5200
HP Laserjet 1320	HP Laserjet Pro MFP M426FDN
HP Laserjet P2015DN	HP Officejet 6000
HP Laserjet P2030 USB Printer	HP Officejet Pro 8000 Enterprise
HP Laserjet P2055	Epson Stylus D120
HP Laserjet Pro MFP M426FDN	Konica Minolta bizhub C3100P
HP Officejet 6000	
HP Photosmart C4680	
HP Photosmart C4780	
Konica Minolta bizhub C3100P	
Lexmark E250D	
Lexmark E460dn	

PCL compatible printers up to PCL version 5 should work properly fine with the M series instruments. Please recheck whether your printer is PCL compatible before you connect it to the instrument.

**If you are using accessories that are not supplied by Anton Paar, we do not guarantee functionality and safety.**