


# AP Connect Editions

DOWNLOAD  
FULL 60-DAY DEMO



[www.anton-paar.com/apc-demo](http://www.anton-paar.com/apc-demo)

	Start (free)	Standard	Pharma
	↓	↓	↓
	For all customers who want to familiarize themselves and grow with AP Connect, and want to connect a single instrument.	Broad-based. For labs and applications across industries.	Focused on labs in intensely regulated environments. For operators with specific stringent testing requirements.
Especially for these industries	<ul style="list-style-type: none"> <li>- Food &amp; beverages</li> <li>- Petroleum</li> <li>- Chemical</li> <li>- Flavors &amp; fragrances</li> </ul>	<ul style="list-style-type: none"> <li>- Food &amp; beverages</li> <li>- Petroleum</li> <li>- Chemical</li> <li>- Flavors &amp; fragrances</li> </ul>	<ul style="list-style-type: none"> <li>- Pharmaceutical</li> <li>- Personal care</li> <li>- Healthcare</li> <li>- Nutraceutical</li> <li>- Government compliance</li> </ul>
License	Free	One-time license fee or annual subscription fee	One-time license fee or annual subscription fee
Instrument connectivity	1 instrument	Up to 30 instruments	Up to 30 instruments
Data collection	✓	✓	✓
Data analysis (results browser)	✓	✓	✓
Data export	✓	✓	✓
Configurable automatic exports	1 rule	Up to 100 rules	Up to 100 rules
Supported file formats for data export	PDF, CSV	PDF, XLS, CSV	PDF, XLS, CSV
Desktop installation	✓	✓	✓
Server installation	✓	✓	✓
Client installation (multi-user access)	1 client	✓	✓
User management supporting Active Directory	✓	✓	✓
Fileless data transfer (REST API) <sup>1</sup>		✓	✓
Connect to instruments from other vendors <sup>2</sup>		✓	✓
Control of instruments (automated measurement series)		✓	✓
Roles and permissions			✓
Review & approval process			✓
Electronic signature			✓
Audit trail			✓
Compliant with Title 21 CFR Part 11			✓

As per the features shown in the table, compliance with the ALCOA++ data integrity principle is ensured.

<sup>1</sup> Via separate license (extra license costs apply).

<sup>2</sup> Additional adapters needed.

© 2023 Anton Paar GmbH | All rights reserved.  
Specifications subject to change without notice.  
D961P008EN-LTR-D



## Software for Your Paperless Lab

AP Connect





# AP Connect.

## Let the Data Flow.

Data can be complex, but this is easy.

Connect and collect, plug and play. All your lab data – just a click away. AP Connect is a lab execution software that connects instruments from Anton Paar and other vendors, communicates measurement information, and ensures compliance. It's your data hub in the lab – your world right at your fingertips.



# Paperless.


## Data at Your Fingertips.

Why is AP Connect different than anything else? The versatile system architecture allows a plug-and-play integration within a few minutes. No complex, time-consuming implementation. No bells and whistles. Just time-saving functionality and increased productivity.

Tens of thousands of measurements in one digital space. 100 % efficiency. Fun, even. Catalyze your lab work.

### Reward your data – in three easy steps

- 1 INSTALL**  
Install AP Connect and, optionally, the instrument adapters in as little as four minutes
- 2 LOG IN**  
Log in to AP Connect in seconds with your Windows access credentials (Active Directory)
- 3 REGISTER**  
Register your instruments in less than one minute



← **Do my instruments work with AP Connect?**  
See all of the 60+ Anton Paar instruments that connect. See all supported instruments from other vendors connected via instrument adapters

[www.anton-paar.com/apc](http://www.anton-paar.com/apc)

Connect.	Collect.	Check.
<ul style="list-style-type: none"><li>→ Plug and play</li><li>→ Connect instruments from Anton Paar and other vendors</li></ul>	<ul style="list-style-type: none"><li>→ One central database for all your measurement data</li><li>→ Access from across the company network</li><li>→ Work with information instead of searching</li></ul>	<ul style="list-style-type: none"><li>→ Call up your data anytime to make faster decisions</li><li>→ Audit trail function for full traceability in your lab work</li><li>→ Configurable approval process of all your data with electronic signature</li></ul>

<b>Zero transcription errors ensure data quality.</b>	<b>Zero paperwork saves time and frees up resources for data analysis.</b>	<b>All data for analysis in one place enables fast and safe decisions.</b>
-------------------------------------------------------	----------------------------------------------------------------------------	----------------------------------------------------------------------------

Control.	Communicate.	Comply.
<ul style="list-style-type: none"><li>→ User management using the Active Directory</li><li>→ Secure assignment of roles &amp; permissions</li><li>→ Trigger measurements with 10+ parameters directly from your PC</li></ul>	<ul style="list-style-type: none"><li>→ One single interface for all your instruments to your data management system</li><li>→ Remote data access: for fileless data transfer via the standardized REST interface</li><li>→ File-based automatic and manual exports</li></ul>	<ul style="list-style-type: none"><li>→ The combination of all AP Connect features ensures full compliance in regulated environments (e.g., FDA 21 CFR Part 11).</li></ul>

<b>Batch measurement increases sample throughput.</b>	<b>Configure your communication paths once, and your data will flow freely.</b>	<b>System audit trail and permission management guarantee full traceability.</b>
-------------------------------------------------------	---------------------------------------------------------------------------------	----------------------------------------------------------------------------------