


Chemical Starter Kit

Buy Both Benchtop Instruments
and Get the Handheld Density
Meter For Free

FIND OUT MORE



DMA 501/DMA 1001 | Abbemat 3001/Abbemat 3101/Abbemat 3201 | DMA 35 Ex

Use Anton Paar's density meter and refractometer to
determine physical properties in the chemical industry.
Precision guaranteed!





1



2



3

Our three-year warranty is available to all customers. All new instruments will include repair for three years.

1

DMA 501/DMA 1001

DMA 501 and DMA 1001 are entry-level digital density meters for the lab and quality checks at the production line and storage facilities. Speed up the measurement of pasty, inhomogeneous, sedimenting, and particle-containing samples.

- 3-digit/4-digit density meters from the market leader
- Filling volume determination of packages by mass to volume conversion
- QC near the production line

2

Abbemat 3001/Abbemat 3101/Abbemat 3201

Budget-friendly benchtop refractometers providing best-build quality, precision, and temperature control.

- Speeding up of binary-solution concentration measurements due to a simplified user interface and predefined methods
- Scratch-proof sapphire prism for long-term precise results
- Easy printing and documentation of results

3

DMA 35 Ex

DMA 35 Ex is the most-sold intrinsically safe portable density meter on the market (Ex marking [Ex] II 2 G Ex ib IIB T4) to measure the density and concentration of your sample directly on-site. A single portable density meter replaces all glass hydrometers and titrators.

- For measuring flammable chemicals in hazardous environments
- Leak-proof sealed housing to withstand rough conditions
- Protection class IP54